






KerraMax Care® Super-absorbent dressing		Product Code	Size	Dressings/Box	HCPCS
 <p>KerraMax Care is a super-absorbent dressing that combines a unique wicking layer to distribute exudate throughout the dressing and absorptive polymers that absorb and retain high levels of exudate into a durable yet conformable dressing. This results in less nursing time spent on dressing changes as well as a reduced risk of maceration.</p>		PRD500-025	2 x 2 in	10	A6196
		PRD500-050	4 x 4 in	10	A6196
		PRD500-100	5 x 6 in	10	A6197
		PRD500-120	4 x 9 in	10	A6197
		PRD500-240	8 x 9 in	10	A6197
		PRD500-380	8 x 12 in	5	A6198
		PRD500-600	8 x 20 in	5	A6198
		PRD500-300	8 x 9 in Multisite	5	A6197

KerraFoam® Gentle Border – Silicone retention dressing					
 <p>KerraFoam® Gentle Border's absorbent core and foam pad absorbs and locks in exudate. The soft silicone adhesive contact layer helps to prevent painful dressing changes and skin stripping making this dressing ideal for both wound care and skin protection.</p>		CWL1010	3 x 3 in	10	A6212
		CWL1011	4 x 4 in	10	A6212
		CWL1013	5 x 5 in	10	A6212
		CWL1132	6 x 6 in	10	A6213
		CWL1135	7 x 8 in Oval	10	A6213
		CWL1012	4 x 8 in	10	A6212
		CWL1133	4 x 10 in	10	A6213
		CWL1134	4 x 12 in	5	A6213
		CWL1023	6.7 x 6.9 in Small Sacrum	5	A6212
		CWL1040	9 x 10 in Large Sacrum	5	A6213
		CWL1024	9.8 x 9 in Heel	5	A6213

KerraFoam® Non-Adhesive					
 <p>KerraFoam® Non-adhesive has the same excellent fluid retention capabilities as KerraFoam® Gentle Border without the silicone wound interface or border.</p>		CWL1116	4 x 5 in	10	A6210
		CWL1117	6 x 6 in	10	A6210
		CWL1118	4 x 8 in	10	A6210
		CWL1119	8 x 8 in	10	A6211
		CWL1131	8 x 12 in	10	A6211

KerraCel® Ag Gelling Fiber dressing – 100% CMC dressing with AG Oxysalts™ Technology					
 <p>KerraCel® Ag is a highly absorbent, conformable silver wound dressing made of 100% CMC containing unique Ag Oxysalts™ Technology. Each Ag Oxysalt ion is up to 6 times more powerful than traditional silver which leads to a faster and more powerful bacteria kill - including superbugs - and the ability to destroy biofilms.</p>		CWL1157	2 x 2 in	10	
		CWL1158	4 x 5 in	10	
		CWL1159	6 x 6 in	5	
		CWL1160	8 x 12 in	5	
		CWL1168	0.75 x 12 in Ribbon	5	
		CWL1162	0.75 x 18 in Ribbon	5	
		CWL1163	16 x 16 in	5	

KerraCel® Gelling Fiber dressing – 100% carboxymethyl cellulose dressing					
 <p>KerraCel® is a highly absorbent, conformable wound dressing made of 100% CMC. It locks in exudate to help protect periwound skin from maceration while contouring to wound bed.</p>		CWL1032	2 x 2 in	10	A6196
		CWL1033	4 x 5 in	10	A6197
		CWL1034	6 x 6 in	5	A6197
		CWL1035	1 x 18 in Ribbon	5	A6199

KerraContact® Ag dressing - with Ag Oxysalts™ Technology



KerraContact® Ag is a contact layer that contains unique Ag Oxysalts™ Technology. Each Ag Oxysalt ion is up to 6 times more powerful than traditional silver which leads to a faster and more powerful bacteria kill - including superbugs - and the ability to destroy biofilms.

Product Code	Size	Dressings/Box	HCPCS
PN-09-0102	2 x 2 in	10	A6206
PN-09-0103	4 x 5 in	10	A6207
PN-09-0104	4 x 8 in	10	A6207
PN-09-0105	6 x 6 in	10	A6207
PN-09-0106	6 x 9 in	10	A6208
PN-09-0108	8 x 16 in	10	A6208
PN-09-0116	16 x 16 in	5	A6208
PN-09-0120	4 x 48 in	5	A6208
PN-09-0132	16 X 32 in	5	A6208

KerraLite Cool® - Hydrogel moisture-balancing dressing



KerraLite Cool® is a soothing, debriding and moisturizing dressing that provides the ideal environment for managing dry to lightly exuding sloughy wounds. It's strong, transparent hydrogel is cool and soothing to the touch and will not dry out or stick to the underlying tissue on removal.

Non-Adhesive			
CWL1004	2.4 x 2.4 in	5	A6234
CWL1005	5 x 3.5 in	5	A6235
CWL1006	7 x 5 in	5	A6235
Adhesive Border			
CWL1007	3 x 3 in	5	A6237
CWL1008	4 x 4 in	5	A6237
CWL1009	6 x 6 in	5	A6238

KerraPro® Pressure reducing pads



Shaped pads are made from 100% silicone and have the ability to redistribute pressure to protect the skin on bony prominences or under medical devices such as masks and collars.

KPRO10	Sheet 4 x 4 x 0.1 in	5	
KPRO20	Sheet 4 x 4 x 0.5 in	5	
KPRO30	Heel/Elbow	2	
KPRO40	Sacrum/Ankle	1	
KPRO50	Strips 12 x 2 x 0.1 in	5	
KPRO60	Strips 20 x 1 x 0.1 in	5	

DebriMitt™ - Gentle, fast wound preparation mitt



DebriMitt™ is a unique polyester-fiber pad used to prepare the wound bed and surrounding skin. It gently removes non-viable tissue, hyperkeratotic skin, and debris in a manner that increases patient tolerance and acceptance with wound preparation.

CWL1021	4 x 4 in	5	
CWL1022	5 x 8 in	5	

Touchless Care® No-Touch skin care products



This patented spray-on line of products allow for easy application to increase efficiency and patient compliance. The Zinc Oxide spray soothes and protects irritated skin with 25% Zinc Oxide and 20% Dimethicone without the mess and pain associated with traditional rubbed-in ointments.

62402	Zinc Oxide Protectant Spray	2.0 oz	24
62404	Zinc Oxide Protectant Spray	4.5 oz	24
72402	Clear Protectant Spray	2.0 oz	24
82402	Antifungal Spray	2.0 oz	24

To place an order please contact your medical supply distributor or call 855-522-2211



Crawford Healthcare Inc. | 2005 South Easton Rd. | Suite 203 | Doylestown, PA 18901
Tel 855-522-2211 | Email us.info@crawfordhealthcare.com
www.crawfordhealthcare.us